

# FOAM PLATFORM SET #REF EPM02


The examination on a **foam platform**, which shunts the podal afferents by blurring the proprioceptive information, directs the postural diagnosis to the cephalic sensors: mandibular, ocular and / or vestibular. It also accentuates pelvic pedal asymmetries.

It can be used for various applications:

- ✔ Anteflexion test / Standing flexion test , the so-called Bassani / Vorlauf test
- ✔ Test of alternating nuchal rotation
- ✔ Osteopathic tests (interspinous in extension, lateral extension at the level of the offending facet, ribs in rotation and inspiration)
- ✔ Visualization of the postural imbalance of the feet and knees
- ✔ Proprioceptive rehabilitation work

The set consists of a platform with an integrated foam and **a set of 2 spare foams**.

# FOAM PLATFORM



## FOAM PLATFORM & RECHARGES

### DESCRIPTION

The integrated foam of high density, has been selected and validated by the most demanding practitioners after more than thirty years of research and clinical use.  
The technical rectangle is adjusted to a protective frame which allows great versatility of use in complete safety!

### FEATURES

- ✓ CARRYING OUT NUMEROUS TESTS
- ✓ LARGE POSITIONING SURFACE
- ✓ OPTIMAL FRAME STABILITY WITH ANTI-SLIP PADS
- ✓ LONG DURABILITY EVEN IN INTENSIVE USE

### TECHNICAL CHARACTERISTICS

dimensions	526 x 436 x 55mm
Weight	4,2kg
Material	EVA50 foam, CDF
Color	Black and white
Maximum load (kg)	130 kg

### CERTIFICATES

- ✓ MADE IN FRANCE
- ✓ 1 YEAR OF USE PER SIDE