# MADDOX TEST SET #REF ETM03

**The Maddox** rod test and the cover test make it possible to orient the diagnosis towards the visual sensor's asynchrony of and also to measure the impact of this sensor's postural treatment.

The set includes a double rod, with an occluder end and a Maddox end, as well as a light point with a ruler.

# **MADDOX TEST SET**



## **DOUBLE ROD: MADDOX + OCCLUDER**

#### **DESCRIPTION**

The Maddox end has a specific red grooved lens that transforms a light point into a line and the other end an occluder of a suitable diameter.

#### **FEATURES**

DIMENSIONS: LENGTH 220 X DIAMETER 60 MM

## LIGHT POINT WITH RULER

#### **DESCRIPTION**

This light source is used specifically with the Maddox rod

The light comes in the form of a round contemporary design box with a magnet, easily adjustable on a metal ruler.

### **FEATURES**

2 IN 1 ROD

PUNCTIFORM VALIDATED LIGHT

SIMPLE TOUCH ON/OFF SWITCH

QUICK HEIGHT ADJUSTMENT OF THE LIGHT

EASY WALL MOUNTING

### TECHNICAL FEATURES

Ruler dimensions	500 mm x 15mm
Weight	0.3 kg
Material	Pvc, steel
Lighting	LED
Supply	3 batteries AAA (not included)
Color	white
Ruler's wall mounting	Adhesive tape + holes